Installation Manual

RadiForce® RX350

Color LCD Monitor

Important

Please read this Installation Manual, and “Instructions for Use” (separate volume) carefully to familiarize yourself with safe and effective usage.

-
- Please refer to the “Instructions for Use” for basic information ranging from connection of the monitor to a PC to using the monitor.
 - The latest Installation Manual is available for download from our web site:
<http://www.eizoglobal.com>
-

This product has been adjusted specifically for use in the region to which it was originally shipped. If operated outside this region, the product may not perform as stated in the specifications.

No part of this manual may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, or otherwise, without the prior written permission of EIZO Corporation. EIZO Corporation is under no obligation to hold any submitted material or information confidential unless prior arrangements are made pursuant to EIZO Corporation's receipt of said information. Although every effort has been made to ensure that this manual provides up-to-date information, please note that EIZO monitor specifications are subject to change without notice.

CONTENTS


CONTENTS	3		
Chapter 1 Basic Operations and Functions	4		
1-1. Basic Adjustment menu operations	4		
1-2. Functions	6		
Chapter 2 Adjusting Screens	7		
2-1. Compatible Resolutions	7		
2-2. Setting the Resolution	7		
2-3. Adjusting the Color.....	8		
● Selecting the display mode (CAL Switch Mode)	8		
● Performing advanced adjustments	9		
● Adjustable settings in each mode	9		
● Adjusting Brightness	9		
● Enabling or disabling Hybrid Gamma PXL	10		
● Adjusting Contrast.....	10		
● Adjusting Temperature	10		
● Adjusting Gamma.....	11		
● Adjusting Hue	11		
● Adjusting Saturation	11		
● Adjusting Gain.....	12		
2-4. Selecting Screen Size	12		
Chapter 3 Setting the Monitor	13		
3-1. Switching among Input Signals	13		
3-2. Setting Input Signal Switching Procedure ...	13		
3-3. Skipping Display Modes That Are Not Used.....	14		
3-4. Setting Language.....	14		
3-5. Showing Monitor Information.....	14		
3-6. Displaying or Hiding the EIZO Logo	15		
3-7. Locking Operation Switches	15		
3-8. Rotating the Image	16		
3-9. Selecting the Input Signal for DDC Communication.....	17		
3-10. Enabling/Disabling DC5V Output.....	17		
3-11. Displaying or Hiding the Grayscale Warning.....	18		
3-12. Enabling/Disabling Sharpness Recovery ...	18		
3-13. Setting the Preferred Refresh Rate.....	19		
3-14. Switching the DisplayPort Version between 1.2 and 1.1.....	19		
3-15. Resetting to Default.....	20		
● Resetting the color adjustments	20		
● Resetting all adjustments to defaults	20		
Chapter 4 Controlling Monitor Quality	21		
4-1. Executing Calibration.....	21		
● Executing SelfCalibration.....	21		
● Viewing SelfCalibration result	21		
4-2. Checking Grayscale.....	22		
● Executing Grayscale Check.....	22		
● Viewing Grayscale Check result	22		
4-3. Acquiring Life Prediction Data	22		
4-4. Setting On/Off for Warning and QC History	23		
Chapter 5 Power Saving Functions	24		
5-1. Setting Power Save.....	24		
5-2. Adjusting the Brightness of Indicator.....	25		
5-3. Setting DisplayPort Power Save	25		
5-4. Setting the Power Saving Mode When a Person Leaves the Monitor.....	26		
Chapter 6 Troubleshooting	27		
6-1. No picture	27		
6-2. Imaging problems	28		
6-3. Other problems	29		
● Error code table.....	30		
Appendix	31		
LIMITED WARRANTY	31		

Chapter 1 Basic Operations and Functions

1-1. Basic Adjustment menu operations

1. Displaying Adjustment Menu

1. Touch one of the switches (excluding ).

2. The guide is displayed.


3. Select  to display the Adjustment menu.


2. Adjusting/Setting

1. Choose a menu to adjust/set with  or , and then choose .

2. Choose an item to adjust/set with **^** or **v** , and then choose **✓**.

3. Adjust/set the selected item with **<** **>** , and then choose **✓**.

3. Exiting

1. Choose **×** a few times. The Adjustment menu finishes.

1-2. Functions

The following table shows all the Adjustment menu's adjustment and setting menus.

Main menu	Item	Reference
CAL Switch Mode	Mode	"2-3. Adjusting the Color" (page 8)
	Brightness	
	Hybrid Gamma PXL	
	Contrast	
	Temperature	
	Gamma	
	Hue	
	Saturation	
	Gain	
Reset		
RadiCS SelfQC	SelfCalibration	"4-1. Executing Calibration" (page 21)
	Grayscale Check	"4-2. Checking Grayscale" (page 22)
	LEA	"4-3. Acquiring Life Prediction Data" (page 22)
	Settings	"4-4. Setting On/Off for Warning and QC History" (page 23)
Signal	Picture Expansion	"2-4. Selecting Screen Size" (page 12)
Preference	Input	"3-1. Switching among Input Signals" (page 13)
	Auto Input Detection	"3-2. Setting Input Signal Switching Procedure" (page 13)
	Power Save	"5-1. Setting Power Save" (page 24)
	Indicator	"5-2. Adjusting the Brightness of Indicator" (page 25)
	Presence Sensor	"5-4. Setting the Power Saving Mode When a Person Leaves the Monitor" (page 26)
	Mode Skip	"3-3. Skipping Display Modes That Are Not Used" (page 14)
	Monitor Reset	"3-14. Switching the DisplayPort Version between 1.2 and 1.1" (page 19)
Languages		"3-4. Setting Language" (page 14)
Information		"3-5. Showing Monitor Information" (page 14)
Administrator Settings	On-Screen Logo	"3-6. Displaying or Hiding the EIZO Logo" (page 15)
	Key Lock	"3-7. Locking Operation Switches" (page 15)
	Image Rotation	"3-8. Rotating the Image" (page 16)
	DDC	"3-9. Selecting the Input Signal for DDC Communication" (page 17)
	DC5V Output	"3-10. Enabling/Disabling DC5V Output" (page 17)
	DP Power Save	"5-3. Setting DisplayPort Power Save" (page 25)
	Grayscale Warning	"3-11. Displaying or Hiding the Grayscale Warning" (page 18)
	Sharpness Recovery	"3-12. Enabling/Disabling Sharpness Recovery" (page 18)
	Signal Format	"3-13. Setting the Preferred Refresh Rate" (page 19) "3-14. Switching the DisplayPort Version between 1.2 and 1.1" (page 19)

*1 The adjustment/settings available on "CAL Switch Mode" depend on the selected mode. (See "Adjustable settings in each mode" (page 9).)

Chapter 2 Adjusting Screens

2-1. Compatible Resolutions

For compatible resolutions of this product, refer to Instruction for Use "Compatible Resolutions".

2-2. Setting the Resolution

Windows 10

1. Right-click the mouse anywhere on the desktop except on icons.
2. From the displayed menu, click "Display Settings".
3. On the "Customize your display" dialog box, click "Advanced display settings".
4. Select a monitor, then choose a resolution from the "Resolution" pull-down menu.
5. Click the "Apply" button.
6. When a confirmation dialog box is displayed, click "Keep changes".

Windows 8.1/Windows 8/Windows 7

1. For Windows 8.1/Windows 8, click the "Desktop" tile from the Start screen.
2. Right-click the mouse anywhere on the desktop except on icons.
3. From the displayed menu, click "Screen resolution".
4. Select a monitor, then choose a resolution from the "Resolution" pull-down menu.
5. Click the "OK" button.
6. When a confirmation dialog box is displayed, click "Keep changes".

Note

- If you want to change the size of displayed characters and other items, select "Display" from Control Panel and change the zoom percentage.
-

Windows Vista

1. Right-click the mouse anywhere on the desktop except on icons.
2. From the displayed menu, click "Personalize".
3. On the "Personalization" window, click "Display Settings".
4. On the "Display Settings" dialog box, select the "Monitor" tab and select desired resolution in the "Resolution" field.
5. Click the "OK" button.
6. When a confirmation dialog box is displayed, click "Yes".

Mac OS X 10.8 and later

1. Select "System Preferences" from the Apple menu.
2. When the "System Preferences" dialog box is displayed, click "Displays". (For Mac OS X 10.8, click "Displays" for "Hardware".)
3. On the displayed dialog box, select the "Display" tab, then choose "Change" for "Resolutions".
4. The list of selectable resolutions is displayed. Select the desired resolution. If the desired resolution is not displayed in the list, hold down the Option key and select "Change".
5. Your selection will be reflected immediately. When you are satisfied with the selected resolution, close the window.

Mac OS X 10.7

1. Select “System Preferences” from the Apple menu.
2. When the “System Preferences” dialog box is displayed, click “Displays” for “Hardware”.
3. On the displayed dialog box, select the “Display” tab, then choose a desired resolution in the “Resolutions” field.
4. Your selection will be reflected immediately. When you are satisfied with the selected resolution, close the window.

2-3. Adjusting the Color

● Selecting the display mode (CAL Switch Mode)

A desired mode can be selected easily according to the monitor application.


Display Modes

√: Calibration is available.

Mode	Purpose	
1-DICOM	Applies gradation characteristics compliant with DICOM Part 14.	√
2-CAL1	Displays the screen adjusted by calibration software.	√
3-CAL2		
4-Custom	Selected for defining a desired set of parameter settings.	-
5-sRGB	Suitable for color matching with sRGB compatible peripherals.	-
6-Text	Suitable for displaying texts for word processing or spreadsheets.	-

Procedure

1. Touch one of the switches (excluding ⏻).
The guide is displayed.
2. Select .
The mode menu appears at the bottom right of the screen.

3. The mode among the list is highlighted in turn each time  is pressed.
You can switch the mode with  or  while the mode menu is displayed.

Note

- The Adjustment menu and the Mode menu cannot be displayed at the same time.

● Performing advanced adjustments

“CAL Switch Mode” of the Adjustment menu allows you to configure and save a separate set of color settings for each mode.

Attention

- As it takes about 15 minutes for the performance of electrical parts to stabilize, adjust the monitor 15 minutes or more after the monitor power has been turned on.
- The same image may be seen in different colors on multiple monitors due to their monitor-specific characteristics. Make fine color adjustment visually when matching colors on multiple monitors.

Note

- The values shown in “%” or “K” are available only as reference.

● Adjustable settings in each mode

The adjustable settings depend on the mode. (You cannot select any non-adjustable or non-settable function.)

√: Adjustable, - : Unadjustable

Function	CAL Switch Mode					
	1-DICOM	2-CAL1	3-CAL2	4-Custom	5-sRGB	6-Text
Brightness	√	√	√	√	√	√
Hybrid Gamma PXL	√	√	√	-	-	-
Contrast	-	-	-	√	√	√
Temperature	-	-	-	√	√	√
Gamma	-	-	-	√	√	√
Hue	-	-	-	√	√	√
Saturation	-	-	-	√	√	√
Gain	-	-	-	√	√	√
Reset	√	√	√	√	√	√

● Adjusting Brightness

The screen brightness is adjusted by changing the brightness of the backlight (light source from the LCD back panel).

Setting Range

0 to 100%

Procedure

1. Select "CAL Switch Mode" of the Adjustment menu, and choose .
2. Select "Brightness" from "CAL Switch Mode", and choose .
3. Use or for adjustment.
4. Select after completing the setting.

Note

- If you think the image is too dark even when the brightness is set to 100%, adjust the contrast.

● Enabling or disabling Hybrid Gamma PXL

When the Hybrid Gamma PXL function is enabled, this product automatically differentiates between monochrome and color parts of the same image at a pixel level, and displays a color image using the gamma value for color image display.

Setting value

On, Off

Procedure

1. Select "CAL Switch Mode" from the Setting menu, then select .
2. Select "Hybrid Gamma PXL", then .
3. Select "On" or "Off", then select .

● Adjusting Contrast

The luminance of the screen is adjusted by varying the video signal level.

Setting Range

0 to 100%

Procedure

1. Select "CAL Switch Mode" of the Adjustment menu, and choose .
2. Select "Contrast" from "CAL Switch Mode", and choose .
"Contrast" appears.
3. Use or for adjustment.
4. Select after completing the setting.

Note

- In the contrast of 100%, every color gradation is displayed.
 - When adjusting the monitor, it is recommended to perform the brightness adjustment which may not lose the gradation characteristics, prior to the contrast adjustment.
 - Perform the contrast adjustment in the following case.
 - If you think the image is too bright even when the brightness is set to 0% (Adjust the contrast to 100% or less.)
-

● Adjusting Temperature

The color temperature can be adjusted.

The color temperature is normally used to express the hue of "white" and/or "black" by a numerical value. The value is expressed in degrees "K" (Kelvin).

The screen becomes reddish at a low color temperature, and bluish at a high color temperature, like the flame temperature. The gain preset values are set for each color temperature setting value.

Setting Range

Native, 6000K to 15000K (in increments of 100K), sRGB, User

Procedure

1. Select "CAL Switch Mode" of the Adjustment menu, and choose .
2. Select "Temperature" from "CAL Switch Mode", and choose .
3. Use or for adjustment.
4. Select after completing the setting.

Note

- "Gain" allows you to perform more advanced adjustments (refer to "Adjusting Gain" (page 12)).
 - If you set "Native", the image is displayed in the preset color of the monitor (Gain: 100% for each RGB).
 - When gain is changed, the color temperature is changed to "User".
-

● Adjusting Gamma

Gamma can be adjusted. The brightness of the monitor varies depending on the input signal, however, the variation rate is not proportional to the input signal. To keep the balance between the input signal and the brightness of the monitor is called as “Gamma correction”.

Setting Range

1.6 to 2.7, sRGB

Procedure

1. Select "CAL Switch Mode" of the Adjustment menu, and choose .
2. Select “Gamma” from “CAL Switch Mode”, and choose .
3. Use or for adjustment.
4. Select after completing the setting.

Attention

- When Hybrid Gamma PXL is set to On, the gamma value for pixels that display a color image is fixed to a value appropriate for color image display regardless of the setting.
-

● Adjusting Hue

The hue can be adjusted.

Setting Range

-100 to 100

Procedure

1. Select "CAL Switch Mode" of the Adjustment menu, and choose .
2. Select “Hue” from “CAL Switch Mode”, and choose .
3. Use or for adjustment.
4. Select after completing the setting.

Attention

- After this adjustment, some color gradations may not be displayed.
-

● Adjusting Saturation

The color saturation can be adjusted.

Setting Range

-100 to 100

Procedure

1. Select "CAL Switch Mode" of the Adjustment menu, and choose .
2. Select “Saturation” from “CAL Switch Mode”, and choose .
3. Use or for adjustment.
4. Select after completing the setting.

Attention

- After this adjustment, some color gradations may not be displayed.
-

Note

- The minimum value (-100) changes the screen to monochrome.
-

● Adjusting Gain

The brightness of each color component red, green, and blue is called Gain. You can change the hue of “white” by adjusting gain.

Setting Range

0 to 100%

Procedure

1. Select "CAL Switch Mode" of the Adjustment menu, and choose .
2. Select "Gain" from "CAL Switch Mode", and choose .
3. Select the color for adjustment from "Red", "Green", and "Blue", and choose .
4. Use or for adjustment.
5. Select after completing the setting.

Attention


- After this adjustment, some color gradations may not be displayed.

Note

- The value of gain changes with that of color temperature.
- When gain is changed, the color temperature is changed to "User".

2-4. Selecting Screen Size

The image with the resolution other than the recommended resolution is displayed in full screen automatically. You can change the screen size by using the "Picture Expansion" function from "Signal". Example: Image size 1024 x 768 (with "Portrait" graphic board setting)

Procedure

1. Select "Signal" of the Adjustment menu, and press .
2. Select "Picture Expansion" from "Signal", and press .
3. Using or , select "Full Screen", "Aspect Ratio", or "Dot by Dot".

Setting	Function
Full Screen	Displays an image in full screen. Images are distorted in some cases because the vertical rate is not equal to the horizontal rate.
Aspect Ratio	Displays an image in full screen. In some cases, a blank horizontal or vertical border appears to equalize the vertical rate and the horizontal rate.
Dot by Dot	Displays images with the specified resolution.

4. Press after completing the setting.


Chapter 3 Setting the Monitor

3-1. Switching among Input Signals

Procedure

1. Select "Preferences" of the Adjustment menu, and choose .
2. Select "Input" from "Preferences", and choose .
3. Using or , select either "DVI" or "DisplayPort".
4. Select after completing the setting.

The selected input port name is displayed.

3-2. Setting Input Signal Switching Procedure

Setting	Function
On	The monitor recognizes the connector through which signals are input. When a PC is turned off or enters power saving mode, the monitor automatically displays another signal.
Off	The monitor displays the signal from the connector currently selected whether or not signals are input. In this case, select the input signal to be displayed by referring to "3-1. Switching among Input Signals" (page 13).

Procedure

1. Select "Preferences" of the Adjustment menu, and choose .
2. Select "Auto Input Detection" from "Preferences", and choose .
3. Select "On" or "Off" with or .
4. Select after completing the setting.

Note

- When the main power supply is turned on or off, signals are automatically detected irrespective of the setting of this function.

3-3. Skipping Display Modes That Are Not Used

When selecting a mode, you can skip particular modes.

This setting is suitable when you are using limited display modes and you do not want to change the customized display state.

Procedure

1. Select "Preferences" of the Adjustment menu, and choose .
2. Select "Mode Skip" from "Preferences", and choose .
3. Select the mode to change with or , and choose .
4. Select "Skip" or "-" with or .
5. Select after completing the setting.

Attention

- You cannot set all the modes disable. Set at least one mode to "-".

3-4. Setting Language

A language for the Adjustment menu and messages can be selected.

Selectable languages

English/Deutsch/Français/Español/Italiano/Svenska/Japanese/Simplified Chinese/Traditional Chinese

Procedure

1. Select "Languages" of the Adjustment menu, and choose .
2. Select a language with or .
3. Select after completing the setting.

3-5. Showing Monitor Information

The information about the current input signals and this product can be displayed.

Procedure

1. Select "Information" of the Adjustment menu, and choose .

The "Information" is displayed.

The input signal information, model name, serial number, firmware version, and usage time can be displayed.

(Example)

Information	
RadiForce RX350	S/N: 12345678
Version	****-****-****
Usage Time (h)	Product: 1234567890
	Backlight: 1234567890
Asset Tag Number	*****
Input Signal	DisplayPort
	1280 X 1024
	fH: 63.79 kHz
	fV: 59.84 Hz
	fD: 108.0 MHz


Attention

- The usage time of the monitor is not always "0" at the time when you purchase it due to the operations such as factory inspection.

3-6. Displaying or Hiding the EIZO Logo

The EIZO logo can be displayed or hidden.

Procedure

1. Touch  to turn off the monitor power.
2. While holding down the switch at the far left () , touch  for at least two seconds to turn on the monitor power.

The "Administrator Settings" menu appears.


3. Select "On-Screen Logo" from "Administrator Settings", and choose .
4. Select "On" or "Off" with  or , and choose .
5. Select "Apply" with  or .
6. Press .

The "Administrator Settings" menu is closed.

3-7. Locking Operation Switches

The adjusted/set status can be locked to prevent changing them.

Procedure

1. Touch  to turn off the monitor power.
2. While holding down the switch at the far left () , touch  for at least two seconds to turn on the monitor power.

The "Administrator Settings" menu appears.

3. Select "Key Lock" from "Administrator Settings", and choose .
4. Using  or , select "Off", "Menu", or "All", and choose .

Setting	Lockable switches
Off (Default setting)	None (All switches are usable)
Menu	Adjustments/settings using the Adjustment menu
All	All switches other than 

5. Select "Apply" with  or .
6. Select .

The "Administrator Settings" menu is closed.

3-8. Rotating the Image

The display direction of the image can be changed according to your environment.

Attention



- Refer to the User's Manual of the graphics board.
- The setting is available even when the PC is not started up.

Procedure

1. Touch  to turn off the monitor power.
2. While holding down the switch at the far left () , touch  for at least two seconds to turn on the monitor power.

The "Administrator Settings" menu appears.

3. Select "Image Rotation" from "Administrator Settings", and choose .
4. Select "Portrait" or "Landscape" with  or , and choose .
5. Select "Apply" with  or .
6. Select after completing the setting.

Monitor orientation	Display example
Landscape Selected for horizontal installation.	
Portrait Selected for vertical installation.	

3-9. Selecting the Input Signal for DDC Communication

Specify the input signal with which DDC communication is used. Use Auto normally. By using DDC communication, you can use RadiCS LE without connecting a USB cable. (This is different from the DDC/CI communication.)

Procedure

1. Touch  to turn off the monitor power.
2. While holding down the switch at the far left () touch  for at least two seconds to turn on the monitor power.

The "Administrator Settings" menu appears.

3. Select "DDC" from "Administrator Settings", and choose .
4. Select "Off", "DVI", "DisplayPort" or "Auto" with  or , and choose .
5. Select "Apply" with  or .
6. Select .

The "Administrator Settings" menu is closed.

Note

- When "Auto" has been selected, the DDC communication of ["3-1. Switching among Input Signals" \(page 13\)](#) is enabled.
-

3-10. Enabling/Disabling DC5V Output

The DC5V output can be enabled/disabled. Turn off this setting normally.

Note

- When using TDL3600 (optional item), turn On this setting.
-

Procedure

1. Select  to turn off the monitor power.
2. While holding down the switch at the far left () touch  for at least two seconds to turn on the monitor power.

The "Administrator Settings" menu appears.

3. Select "DC5V Output" from "Administrator Settings", and choose .
4. Select "On" or "Off" with  or , and choose .
5. Select "Apply" with  or .
6. Select .

The "Administrator Settings" menu is closed.


3-11. Displaying or Hiding the Grayscale Warning

When the DisplayPort signal is input in a recommended resolution, the grayscale display error message may be displayed if the monitor detects a 6-bit (64-gradation) signal. When the error message is displayed, turn off the main power of the monitor, and then turn it on again.

Attention

- Turn "On" this setting normally (the factory default is "On").
 - Turn "Off" this setting only when you use the monitor in 6-bit (64-gradation).
-


Procedure

1. Select  to turn off the monitor power.
 2. While holding down the switch at the far left () , touch  for at least two seconds to turn on the monitor power.
The "Administrator Settings" menu appears.
 3. Select "Grayscale Warning" from "Administrator Settings", and choose .
 4. Select "On" or "Off" with  or , and choose .
 5. Select "Apply" with  or .
 6. Select .
- The "Administrator Settings" menu is closed.

3-12. Enabling/Disabling Sharpness Recovery

Sharpening image can be enabled or disabled. Enabling it enhances the sharpness of images.

Procedure

1. Select  to turn off the monitor power.
 2. While holding down the switch at the far left () , touch  for at least two seconds to turn on the monitor power.
The "Administrator Settings" menu appears.
 3. Select "Sharpness Recovery" from "Administrator Settings", and choose .
 4. Select "On" or "Off" with  or , and choose .
 5. Select "Apply" with  or .
 6. Select .
- The "Administrator Settings" menu is closed.


Attention

- After changing the setting, perform a visual checkup based on medical standards/guidelines.
-


3-13. Setting the Preferred Refresh Rate

The priority of refresh rates can be defined. The "Low" refresh rate is set to 46 Hz (or 47 Hz for landscape), and "High" is to 60 Hz.

Procedure

1. Select  to turn off the monitor power.
2. While holding down the switch at the far left () , touch  for at least two seconds to turn on the monitor power.

The "Administrator Settings" menu appears.

3. Select "Signal Format" from "Administrator Settings", and choose .
4. Select "DVI" with  or , and choose .
5. Select "Low" or "High" with  or , and choose .
6. Select "Apply" with  or .
7. Select .

The "Administrator Settings" menu is closed.


3-14. Switching the DisplayPort Version between 1.2 and 1.1

The DisplayPort version can be switched.


Attention

- Use version 1.2 when setting up a daisy-chain connection.
- When a PC is connected to the DisplayPort connector, the screen may not be displayed depending on the PC due to difference of the DisplayPort version. If that occurs, change this setting.

Procedure

1. Select  to turn off the monitor power.
2. While holding down the switch at the far left () , touch  for at least two seconds to turn on the monitor power.

The "Administrator Settings" menu appears.

3. Select "Signal Format" from "Administrator Settings", and choose .
4. Select "DisplayPort" with  or , and choose .
5. Select "1.1" or "1.2" with  or , and choose .
6. Select "Apply" with  or .
7. Select .

The "Administrator Settings" menu is closed.

3-15. Resetting to Default

There are two different reset functions: One that only resets the color adjustments to default settings, and the other that resets all settings to defaults.

Attention

- After resetting, you cannot undo the operation.

Note

- For default settings, refer to Instruction for Use "Main Default Settings".
-

● Resetting the color adjustments

Only the color adjustment value for the currently selected mode can be reset to the default settings.

Procedure

1. Select "CAL Switch Mode" of the Adjustment menu, and choose .
2. Select "Reset" from "CAL Switch Mode", and choose .
3. Select "OK" with or .
4. Select .

The color adjustments revert to the default settings.

● Resetting all adjustments to defaults

All settings are reset to defaults (except the "Input" and "Administrator Settings" menus).

Procedure

1. Select "Preferences" of the Adjustment menu, and choose .
2. Select "Monitor Reset" from "Preferences", and choose .
3. Select "OK" with or .
4. Select .

All settings excluding the "Input" and "Administrator Settings" menus are reset to defaults.

Chapter 4 Controlling Monitor Quality

This product has the built-in Integrated Front Sensor, which allows the user to perform periodic calibration and Grayscale Check for the monitor independently.

Note

- To perform advanced quality control that conforms to medical standards, use the optional monitor quality control software kit “RadiCS UX1”.
- Perform correlation using RadiCS/RadiCS LE to adjust the measurement result of the Integrated Front Sensor to the measurement result of the UX1 sensor for reference. Periodical correlation enables the measurement accuracy of the Integrated Front Sensor equivalent to that of the UX1 sensor. For details, refer to RadiCS/RadiCS LE User’s Manual.
- “RadiCS” allows you to set the adjustment target for the calibration and judgment value for Grayscale Check, and also the schedule for the calibration and Grayscale Check.
- The supplied “RadiCS LE” software also allows you to set the adjustment target and schedule for the calibration. For details, refer to RadiCS LE User’s Manual (on the CD-ROM).

4-1. Executing Calibration

This function allows you to execute SelfCalibration and view the latest result. In addition, when the schedule is set using RadiCS/RadiCS LE, the next calibration schedule appears.

● Executing SelfCalibration

Procedure

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose .
2. Choose “SelfCalibration” from the “RadiCS SelfQC” menu, and choose .
3. Select “Execute” with or .
4. Choose to exit.
SelfCalibration is executed.

Attention

- When the power of the PC is turned off or input signals are switched during execution of SelfCalibration, the execution is canceled.
- After SelfCalibration has been executed, perform the Grayscale Check based on the actual temperature and illuminance of the working environment.

● Viewing SelfCalibration result

Attention

- The result can be viewed when “QC History” is “On”. (See [“4-4. Setting On/Off for Warning and QC History” \(page 23\)](#).)

Procedure

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose .
2. Choose “SelfCalibration” from the “RadiCS SelfQC” menu, and choose .
3. Select “Result” with or .
4. Select “DICOM”, “CAL1” or “CAL2” with or .
5. Choose to exit.
The SelfCalibration result (Max. Error Rate, Last Time, and Error Code) appears.

4-2. Checking Grayscale

This function allows you to execute Grayscale Check and view the latest result. In addition, when the schedule is set using RadiCS, the next grayscale check schedule appears.

● Executing Grayscale Check

Attention

- Perform the Grayscale Check based on the actual temperature and illuminance of the working environment.
-

Procedure

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose .
2. Choose “Grayscale Check” from the “RadiCS SelfQC” menu, and choose .
3. Select “Execute” with or .
4. Choose to exit.
Grayscale Check is executed.

● Viewing Grayscale Check result

Attention

- The result can be viewed when “QC History” is “On”.
-

Procedure

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose .
2. Choose “Grayscale Check” from the “RadiCS SelfQC” menu, and choose .
3. Select “Result” with or .
4. Select “DICOM”, “CAL1” or “CAL2” with or .
5. Choose to exit.
Up to the five latest Grayscale Check results (Passed/Failed/Canceled/Error) appear.

4-3. Acquiring Life Prediction Data

The data necessary to predict the monitor life can be acquired.

Note

- Data are collected every 100 hours.
 - When you select “Power Save”, data are collected when the mode switches to Power Save or when the power is shut down after 100 hours.
 - When you select “Routine” or “Power Save”, the Integrated Front Sensor performs measurement during data collection. When you select “Off”, measurement using the Integrated Front Sensor is not performed during data collection.
 - Data are saved after the following times have elapsed.
 - 500, 1000, 2000, 4000, 7000, 10000, 15000, 20000, 25000, 30000 hours
-

Procedure

1. Select “RadiCS SelfQC” from the Adjustment menu, and choose .
2. Select “LEA” from “RadiCS SelfQC”, and choose .
3. Select “Meas. Frequency” with or .
4. Select “Power Save”, “Routine”, or “Off” with or .
5. Choose to exit.

4-4. Setting On/Off for Warning and QC History

This function allows you to set On/Off for Warning and QC History when the Grayscale Check result is Failed.

Procedure

1. Choose "RadiCS SelfQC" from the Adjustment menu, and choose .
2. Choose "Settings" from the "RadiCS SelfQC" menu, and choose .
3. Select "Warning" or "QC History" with or .
4. Select "On" or "Off" with or .
5. Choose to exit.

Chapter 5 Power Saving Functions

5-1. Setting Power Save

The monitor can be set to the power saving mode according to the PC status. When the monitor has shifted to the power saving mode, the images on the screen are not displayed.

Attention

- When not using the monitor, you can turn off the main power supply or disconnect the power plug so that the power is cut completely.
- Devices connected to the USB upstream port and USB downstream port work when the monitor is in power saving mode or when the power button of the monitor is Off. Therefore, power consumption of the monitor varies with connected devices even in the power saving mode.
- Five seconds before the monitor enters the power saving mode, it displays a preliminary message.

Procedure

1. Select “Preferences” of the Adjustment menu, and choose .
2. Select “Power Save” from “Preferences”, and choose .
3. Select “Off”, “High”, or “Low” with or .

Setting	Function
Low	Turns down the backlight to the lowest setting. The time until the operation stabilizes after exiting the Power Save mode next time can be reduced.
High	Turns off the backlight. The power saving effect is the greatest.

4. Select after completing the setting.

Power Saving System

At DVI signal input

This monitor complies with the “DVI DMPM” standard.

At DisplayPort signal input

This monitor complies with the “DisplayPort Standard V1.2a”.

The monitor enters the power saving mode after five seconds in connection with the PC setting.

PC	Monitor	Power indicator
On	Operation mode	Green
Power saving mode	Power saving mode	Orange

5-2. Adjusting the Brightness of Indicator

The brightness of the power indicator (green) when the screen is displayed can be set. (By default, the power indicator lights up upon power-on and is set to 4.)

Procedure

1. Select "Preferences" of the Adjustment menu, and choose .
2. Select "Indicator" from "Preferences", and choose .
3. Using or , select the indicator brightness "Off" or in the range from 1 to 7 as desired.
4. Select after completing the setting.

5-3. Setting DisplayPort Power Save

If the PC is connected to the DisplayPort connector, and when the power is turned on or off or when the monitor recovers from power saving mode, window or icon positions may be shifted. If that occurs, set this function to "Off".

Procedure

1. Select to turn off the monitor power.
2. While holding down the switch at the far left () , touch for at least two seconds to turn on the monitor power.
The "Administrator Settings" menu appears.
3. Select "DP Power Save" from "Administrator Settings", and choose .
4. Select "On" or "Off" with or , and choose .
5. Select after completing the setting.

5-4. Setting the Power Saving Mode When a Person Leaves the Monitor

The sensor on the front side of the monitor detects the movement of a person by setting Presence Sensor to On. When a person moves away from the monitor, it shifts to power saving mode automatically and does not display the images on the screen to reduce the power consumption. When a person comes near the monitor again, it recovers from the power saving mode and displays the images. The sensitivity, and time until the power saving mode is activated can be set according to the monitor usage environment and movement of the user.

Note

- At the time of shifting to power saving mode, the message that indicates shifting is displayed.

Attention

- This setting cannot be used during “RadiCS SelfQC” execution.

Procedure

1. Select “Preferences” of the Adjustment menu, and choose .
2. Select “Presence Sensor” from “Preferences”, and choose .
3. Select “On” or “Off” with or .
4. Select after completing the setting.
5. If you have selected “On”, configure the detail settings.

Item	Setting Range	Description
Time	5, 30 sec 1, 3, 5, 10, 15, 30, 45, 60 min	Specifies the time from when a person leaves the monitor until the “No presence detected.” message appears. The monitor shifts to power saving mode about 20 seconds after the message is displayed.
Sensitivity	Level 5 Level 4 Level 3 Level 2 Level 1	When the setting level is high (“Level 5” at maximum), even minor movements of a person are detected, causing the monitor to shift to the power saving mode only infrequently. On the other hand, when the setting level is low (“Level 1” at minimum), only major movements of a person are detected, allowing the monitor to shift to power saving mode sooner. The presence sensor on the front side of the monitor functions by detecting the movement of heat. Consequently, the sensitivity is affected by clothes and the ambient temperature. If the monitor enters power saving mode or recovers from it unintentionally, change the settings.
Reset	-	All adjustments are reset to the default settings.



Note

- The sensitivity becomes lower when the ambient temperature is high, and higher when the temperature is low. Accordingly, you can set the sensitivity to a higher level when temperatures are higher, and to a lower level when temperatures are lower to ensure more appropriate operation.

Chapter 6 Troubleshooting

If a problem still remains after applying the suggested remedies, contact your local EIZO representative.

6-1. No picture

Problem	Possible cause and remedy
<p>1. No picture</p> <ul style="list-style-type: none"> Power indicator does not light. <hr/> <ul style="list-style-type: none"> Power indicator is lighting green. <hr/> <ul style="list-style-type: none"> Power indicator is lighting orange. <hr/> <ul style="list-style-type: none"> Power indicator flashes orange and green. 	<ul style="list-style-type: none"> Check whether the power cord is connected properly. Turn the main power switch on. Touch . Turn off the main power, and then turn it on again. <hr/> <ul style="list-style-type: none"> Increase "Brightness", "Contrast", and/or "Gain" in the Adjustment menu (see "Performing advanced adjustments" (page 9)). Turn off the main power, and then turn it on again. <hr/> <ul style="list-style-type: none"> Switch the input signal (see "3-1. Switching among Input Signals" (page 13)). Move the mouse or press any key on the keyboard. Check whether the PC is turned on. When Presence Sensor is set to "On", the monitor may have shifted to power saving mode. Come near to the monitor. Is the signal cable connected to  for inputting the DisplayPort signal.  is used for output when a daisy-chain connection is set up. Turn off the main power, and then turn it on again. <hr/> <ul style="list-style-type: none"> When connected to DisplayPort, switch the version of DisplayPort to 1.1.(see "3-14. Switching the DisplayPort Version between 1.2 and 1.1" (page 19)). The device that is connected to DVI or DisplayPort have a problem. Solve the problem, turn off the monitor, and then turn it on again. Refer to the manual of the output device for details.
<p>2. The message below appears.</p> <ul style="list-style-type: none"> This message appears when no signal is input. Example: <div data-bbox="284 1473 657 1615" data-label="Image"> </div> <hr/> <ul style="list-style-type: none"> The message shows that the input signal is out of the specified frequency range. (Such signal frequency is displayed in magenta.) Example: <div data-bbox="284 1928 657 2069" data-label="Image"> </div>	<p>This message appears when the signal is not input correctly even when the monitor functions properly.</p> <hr/> <ul style="list-style-type: none"> The message shown left may appear, because some PCs do not output the signal soon after power-on. Check whether the PC is turned on. Check whether the signal cable is connected properly. Switch the input signal (see "3-1. Switching among Input Signals" (page 13)). When connected to DisplayPort, switch the version of DisplayPort to 1.1.(see "3-14. Switching the DisplayPort Version between 1.2 and 1.1" (page 19)). Is the signal cable connected to  for inputting the DisplayPort signal.  is used for output when a daisy-chain connection is set up. Turn off the main power, and then turn it on again. <hr/> <ul style="list-style-type: none"> Check whether the PC is configured to meet the resolution and vertical scan frequency requirements of the monitor (see "Compatible Resolutions" in the "Instructions for Use"). Reboot the PC. Select an appropriate setting using the graphics board's utility. Refer to the manual of the graphics board for details.

6-2. Imaging problems

Problem	Possible cause and remedy
1. The screen is too bright or too dark.	<ul style="list-style-type: none">• Use “Brightness” or “Contrast” in the Adjustment menu to adjust it. (The LCD monitor backlight has a limited life span. When the screen becomes dark or begins to flicker, contact your local EIZO representative.)
2. Characters are blurred.	<ul style="list-style-type: none">• Check whether the PC is configured to meet the resolution and vertical scan frequency requirements of the monitor (see “Compatible Resolutions” in the “Instructions for Use”).
3. Afterimages appear.	<ul style="list-style-type: none">• Afterimages are particular to LCD monitors. Avoid displaying the same image for a long time.• Use the screen saver or power saving function to avoid displaying the same image for extended periods of time.
4. Green/red/blue/white dots or defective dots remain on the screen.	<ul style="list-style-type: none">• This is due to LCD panel characteristics and is not a failure.
5. Interference patterns or pressure marks remain on the screen.	<ul style="list-style-type: none">• Leave the monitor with a white or black screen. The symptom may disappear.

6-3. Other problems

Problem	Possible cause and remedy
1. The adjustment menu does not appear.	<ul style="list-style-type: none"> • Check whether the operation lock function works (see “3-7. Locking Operation Switches” (page 15)).
2. No mode menu is displayed.	<ul style="list-style-type: none"> • Check whether the operation lock function works (see “3-7. Locking Operation Switches” (page 15)).
3. The operation switches do not work.	<ul style="list-style-type: none"> • Check whether the operation lock function works (see “3-7. Locking Operation Switches” (page 15)). • Check whether there are any water droplets or foreign objects on the switch surface. Gently wipe the switch surface and try operating the switches again with dry hands. • Are you wearing gloves? If so, remove the gloves and try operating the switches again with dry hands.
4. (When using Presence Sensor) The images remain to be displayed when you move away from the monitor. / The images do not reappear when you come near it.	<ul style="list-style-type: none"> • Check the setting environment of the monitor. Presence Sensor may not work correctly in the following environments. <ul style="list-style-type: none"> - The monitor is located in a place affected by wind. - There is some heat emitting device near the monitor. - There is some obstacle in front of the sensor. • Check whether any dirt on the sensor. Clean the sensor with a soft cloth. • Check whether you are in front of the monitor. The Presence Sensor function is activated when the sensor on the front of the monitor detects the movement of heat. • In the high temperature conditions, the monitor may not recover from the power saving mode. If the screen is not displayed even after moving the mouse or pressing any key on the keyboard, turn the monitor off and on using its power button. When the screen is displayed, increase the detection sensitivity value (see “5-4. Setting the Power Saving Mode When a Person Leaves the Monitor” (page 26)).
5. The monitor connected with the USB cable is not detected. / USB devices connected to the monitor do not work.	<ul style="list-style-type: none"> • Check whether the USB cable is connected correctly (see “Making Use of USB (Universal Serial Bus)” in the “Instructions for Use”). • Change the USB port to another one. If the PC or peripheral devices works correctly by changing the USB port, contact your local EIZO representative. (Refer to the manual of the PC for details.) • Reboot the PC. • If the peripheral devices work correctly when the PC and peripheral devices are connected directly, please contact your local EIZO representative. • Check whether the PC and OS are USB compliant. (For USB compliance of the respective devices, consult their manufacturers.) • Check the PC’s BIOS setting for USB when using Windows. (Refer to the manual of the PC for details.)

Problem	Possible cause and remedy
6. The power turns off during operation and no image is displayed.	<ul style="list-style-type: none"> The monitor power is turned off automatically when the internal temperature is extremely high. Change the installation location or environment, and then turn it on again a few minutes later. <ul style="list-style-type: none"> The ventilation slot must not be closed. There is no high temperature object around. If the power is not turned on even after the installation location or environment has changed, contact your local EIZO representative. <p>Attention</p> <ul style="list-style-type: none"> Do not block the ventilation hole of the monitor. Do not use the monitor in a place with poor ventilation.
7. Integrated Front Sensor stays out.	<ul style="list-style-type: none"> Turn the main power off and on again.
8. It takes time to start RadiCS.	<ul style="list-style-type: none"> Connect the attached USB cable.
9. SelfCalibration/Grayscale Check failure	<ul style="list-style-type: none"> Refer to the error code table in the “Instructions for Use”. If an error code that does not appear in the error code table is displayed, contact your local EIZO representative. <p>Attention</p> <ul style="list-style-type: none"> Do not touch the Integrated Front Sensor.
10. Time is not displayed correctly in the history of SelfCalibration/Grayscale Check.	<ul style="list-style-type: none"> Execute Monitor Detection using RadiCS/RadiCS LE. For details, refer to RadiCS/RadiCS LE User’s Manual.

● Error code table

Error Code	Description
0***	<ul style="list-style-type: none"> Errors that occurred during SelfCalibration
1***	<ul style="list-style-type: none"> Errors that occurred during Grayscale Check
*1**	<ul style="list-style-type: none"> Errors that occurred during DICOM
*2**	<ul style="list-style-type: none"> Errors that occurred during CAL1
*3**	<ul style="list-style-type: none"> Errors that occurred during CAL2
**10	<ul style="list-style-type: none"> The product maximum brightness may be lower than the target brightness. Lower the target brightness.
**11	<ul style="list-style-type: none"> The product minimum brightness may be higher than the target brightness. Raise the target brightness.
**34	<ul style="list-style-type: none"> The sensor may not have come out during calibration, or light may have entered the sensor. Turn off the main power, wait for a few minutes before turning the power back on, and then execute SelfCalibration/Grayscale Check again.
**61	<ul style="list-style-type: none"> The sensor may not have come out. Check whether there is any foreign object near the sensor. Execute SelfCalibration/Grayscale Check again.
**95	<ul style="list-style-type: none"> The execution conditions for SelfCalibration/Grayscale Check may not be appropriate. Check the execution conditions and reset them if necessary.

Appendix

LIMITED WARRANTY

LIMITED WARRANTY

EIZO Corporation (hereinafter referred to as “EIZO”) and distributors authorized by EIZO (hereinafter referred to as the “Distributors”) warrant, subject to and in accordance with the terms of this limited warranty (hereinafter referred to as the “Warranty”), to the original purchaser (hereinafter referred to as the “Original Purchaser”) who purchased the product specified in this document (hereinafter referred to as the “Product”) from EIZO or Distributors, that EIZO and Distributors shall, at their sole discretion, either repair or replace the Product at no charge if the Original Purchaser becomes aware within the Warranty Period (defined below) that (i) the Product malfunctions or is damaged in the course of normal use of the Product in accordance with the description in the instruction manual attached to the Product (hereinafter referred to as the “User’s Manual”), or (ii) the LCD panel and brightness of the Product cannot maintain the recommended brightness specified in the User’s Manual in the course of normal use of the Product in accordance with the description of the User’s Manual.

The period of this Warranty is five (5) years from the date of purchase of the Product (hereinafter referred to as the “Warranty Period”), provided that the brightness of the Product shall be warranted only if the Product has been used within the recommended brightness described in the User’s Manual and the Warranty Period for the brightness is limited to five (5) years from the date of purchase of the Product subject to the usage time being less than or equal to 20,000 hours (the brightness is 500cd/m² and color temperature is 7500K) or less than or equal to 30,000 hours (the brightness is 400cd/m² and color temperature is 7500K). EIZO and Distributors shall bear no liability or obligation with regard to the Product in relation to the Original Purchaser or any third parties other than as provided under this Warranty.

EIZO and Distributors will cease to hold or store any parts (excluding design parts) of the Product upon expiration of seven (7) years after the production of such parts is discontinued. In repairing the monitor, EIZO and Distributors will use renewal parts which comply with our QC standards. If the unit cannot be repaired due to its condition or the stockout of a relevant part, EIZO and Distributors may offer the replacement by a product with equivalent performance instead of repairing it.

The Warranty is valid only in the countries or territories where the Distributors are located. The Warranty does not restrict any legal rights of the Original Purchaser.

Notwithstanding any other provision of this Warranty, EIZO and Distributors shall have no obligation under this Warranty whatsoever in any of the cases as set forth below:

- (a) Any defect of the Product caused by freight damage, modification, alteration, abuse, misuse, accident, incorrect installation, disaster, faulty maintenance and/or improper repair by third party other than EIZO and Distributors;
- (b) Any incompatibility of the Product due to possible technical innovations and/or regulations;
- (c) Any deterioration of the sensor;
- (d) Any deterioration of display performance caused by the deterioration of expendable parts such as the LCD panel and/or backlight, etc. (e.g. changes in brightness uniformity, changes in color, changes in color uniformity, defects in pixels including burnt pixels, etc.);
- (e) Any deterioration of the Product caused by the use at higher brightness than the recommended brightness described in the User’s Manual;
- (f) Any defect of the Product caused by external equipment;
- (g) Any defect of the Product on which the original serial number has been altered or removed;
- (h) Any normal deterioration of the product, particularly that of consumables, accessories, and/or attachments (e.g. buttons, rotating parts, cables, User’s Manual, etc.); and
- (i) Any deformation, discoloration, and/or warp of the exterior of the product including that of the surface of the LCD panel.

To obtain service under the Warranty, the Original Purchaser must deliver the Product, freight prepaid, in its original package or other adequate package affording an equal degree of protection, assuming the risk of damage and/or loss in transit, to the local Distributor. The Original Purchaser must present proof of purchase of the Product and the date of such purchase when requesting services under the Warranty.

The Warranty Period for any replaced and/or repaired product under this Warranty shall expire at the end of the original Warranty Period.

EIZO OR DISTRIBUTORS ARE NOT RESPONSIBLE FOR ANY DAMAGE TO, OR LOSS OF, DATA OR OTHER INFORMATION STORED IN ANY MEDIA OR ANY PART OF ANY PRODUCT RETURNED TO EIZO OR DISTRIBUTORS FOR REPAIR.

EIZO AND DISTRIBUTORS MAKE NO FURTHER WARRANTIES, EXPRESSED OR IMPLIED, WITH RESPECT TO THE PRODUCT AND ITS QUALITY, PERFORMANCE, MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR USE. IN NO EVENT SHALL EIZO OR DISTRIBUTORS BE LIABLE FOR ANY INCIDENTAL, INDIRECT, SPECIAL, CONSEQUENTIAL OR OTHER DAMAGE WHATSOEVER (INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF PROFIT, BUSINESS INTERRUPTION, LOSS OF BUSINESS INFORMATION, OR ANY OTHER PECUNIARY LOSS) ARISING OUT OF THE USE OR INABILITY TO USE THE PRODUCT OR IN ANY CONNECTION WITH THE PRODUCT, WHETHER BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE, EVEN IF EIZO OR DISTRIBUTORS HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. THIS EXCLUSION ALSO INCLUDES ANY LIABILITY WHICH MAY ARISE OUT OF THIRD PARTY CLAIMS AGAINST THE ORIGINAL PURCHASER. THE ESSENCE OF THIS PROVISION IS TO LIMIT THE POTENTIAL LIABILITY OF EIZO AND DISTRIBUTORS ARISING OUT OF THIS LIMITED WARRANTY AND/OR SALES.

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