User's Manual

RadiForce® GS521-BL GS521-CL GS521-CLG

Monochrome LCD Monitor

Important

Please read PRECAUTIONS, this User's Manual, and the Setup Manual (separate volume) carefully to familiarize yourself with safe and effective usage.

-
- Please refer to the Setup Manual for basic information ranging from connection of the monitor to a PC to using the monitor.
 - The latest User's Manual is available for download from our web site:
<http://www.eizo.com>
-

This product has been adjusted specifically for use in the region to which it was originally shipped. If operated outside this region, the product may not perform as stated in the specifications.

- It shall be assured that the final system is in compliance to IEC60601-1-1 requirement.
- Power supplied equipment can emit electromagnetic waves, that could influence, limit or result in malfunction of the monitor. Install the equipment in a controlled environment, where such effects are avoided.
- This is a monitor intended for use in a medical imaging system.

No part of this manual may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, or otherwise, without the prior written permission of EIZO NANA O CORPORATION.

EIZO NANA O CORPORATION is under no obligation to hold any submitted material or information confidential unless prior arrangements are made pursuant to EIZO NANA O CORPORATION's receipt of said information. Although every effort has been made to ensure that this manual provides up-to-date information, please note that EIZO monitor specifications are subject to change without notice.

Apple, Mac, Macintosh, iMac, eMac, Mac OS, MacBook, PowerBook, ColorSync, QuickTime, and iBook are registered trademarks of Apple Inc.

Windows, Windows Vista, Windows Media, SQL Server, and Xbox 360 are registered trademarks of Microsoft Corporation in the United States and other countries.

A round gothic bold bit map font designed by Ricoh is used for the characters displayed on this product.

VESA is a registered trademark or a trademark of the Video Electronics Standards Association in the United States and other countries.

EIZO, the EIZO Logo, ColorEdge, DuraVision, FlexScan, FORIS, RadiForce, RadiCS, RadiNET, Raptor, and ScreenManager are registered trademarks of EIZO NANA O CORPORATION in Japan and other countries.

All other company and product names are trademarks or registered trademarks of their respective owners.

Notice for this monitor

This product is suited to displaying medical images of such modalities as MRI, CT and digital mammography.

This product has been adjusted specifically for use in the region to which it was originally shipped. If the product is used outside the region, it may not operate as specified in the specifications.

This product may not be covered by warranty for uses other than those described in this manual.

The specifications noted in this manual are only applicable when the following are used:

- Power cords provided with the product
 - Signal cables specified by us
-

Only use optional products manufactured or specified by us with this product.

As it takes about 30 minutes for the performance of electrical parts to stabilize, adjust the monitor 30 minutes or more after the monitor power has been turned on.

Monitors should be set to a lower brightness to reduce changes in luminosity caused by long-term use and maintain a stable display.

When the screen image is changed after displaying the same image for extended periods of time, an afterimage may appear. Use the screen saver or power save function to avoid displaying the same image for extended periods of time.

Periodic cleaning is recommended to keep the monitor looking new and to prolong its operation lifetime. (Refer to “[Cleaning](#)” on the next page.)

The screen may have defective pixels or a small number of light dots on the screen. This is due to the characteristics of the panel itself, and is not a malfunction of the product.

The backlight of the LCD panel has a fixed lifetime. When the screen becomes dark or begins to flicker, please contact your dealer.

Do not press on the panel or edge of the frame strongly, as this may result in display malfunctions, such as interference patterns, etc. If pressure is continually applied to the panel, it may deteriorate or damage your panel. (If the pressure marks remain on the panel, leave the monitor with a black or white screen. The symptom may disappear.)

Do not scratch or press on the panel with any sharp objects, as this may result in damage to the panel. Do not attempt to brush with tissues as this may scratch the panel.

When the monitor is cold and brought into a room or the room temperature goes up quickly, dew condensation may occur on the interior and exterior surfaces of the monitor. In that case, do not turn the monitor on. Instead wait until the dew condensation disappears, otherwise it may cause some damage to the monitor.

Cleaning

Attention

- Chemicals such as alcohol and antiseptic solution may cause gloss variation, tarnishing, and fading of the cabinet or panel, and also quality deterioration of the image.
- Never use any thinner, benzene, wax, and abrasive cleaner, which may damage the cabinet or panel.

Note

- The optional ScreenCleaner is recommended for cleaning the cabinet and panel surface. (“ScreenCleaner” is supplied with GS521-CLG.)

If necessary, the stains on the cabinet and panel surface can be removed by moistening part of a soft cloth with water.

To use the monitor comfortably

- An excessively dark or bright screen may affect your eyes. Adjust the brightness of the monitor according to the environmental conditions.
- Staring at the monitor for a long time tires your eyes. Take a 10-minute rest every hour.

CONTENTS

Cover	1	Chapter 4 Power Saving Function	17
Notice for this monitor	3	4-1. Setting the Power Saving [Power Save]	17
CONTENTS	5	4-2. Setting Power Indicator [Power Indicator] ..	17
Chapter 1 Introduction	6	Chapter 5 Troubleshooting	18
1-1. Features.....	6	Chapter 6 Reference	20
1-2. Switches and Indicators.....	6	6-1. Attaching an Arm	20
1-3. Utility Disk.....	7	6-2. Connecting Two PCs to the Monitor	21
● Disk contents and software overview	7	● To switch the input signal.....	21
● To RadiCS LE or ScreenManager Pro for Medical	7	● To set input signal selection [Input Selection]..	22
1-4. Basic Operation and Functions.....	8	6-3. Making Use of USB (Universal Serial Bus). 23	
Basic operation of Adjustment menu.....	8	● Required System Environment	23
Functions	9	● Connection Procedure (Setup of USB Function).	23
Chapter 2 Adjusting Screen	10	6-4. Displaying Monitor Information	24
2-1. Setting Screen Resolution.....	10	● Displaying signal information [Signal Info].....	24
Compatible Resolutions/Frequencies.....	10	● Displaying monitor information [Monitor Info] ..	24
Setting Resolution	10	6-5. Specifications.....	25
● Windows 7.....	10	6-6. Glossary	29
● Windows Vista	10	FCC Declaration of Conformity	30
● Windows XP.....	10	EMC Information	31
● Mac OS X.....	10	LIMITED WARRANTY	34
2-2. Selecting Display Mode.....	11		
2-3. Adjusting Brightness [Brightness].....	11		
Chapter 3 Setting Monitor	12		
3-1. Setting Orientation [Orientation]	12		
3-2. Setting the Input Signals [Signal Selection]13			
● DVI Connection (Dual Link/Single Link)	13		
● DVI Connection (SDG/Packed Pixel).....	13		
● DisplayPort Connection	14		
3-3. Enabling/Disabling Mode Selection [Mode Preset]	14		
3-4. Locking Switches [Key Lock]	14		
3-5. Setting the EIZO Logo Display [Logo]	15		
3-6. Setting Language [Language]	15		
3-7. Setting the Display Position of the Adjustment Menu [Menu Position].....	15		
3-8. Restoring the Default Setting	16		
● To reset brightness adjustment values [Mode Reset].....	16		
● To reset all adjustments to the factory default settings [All Reset]	16		


Chapter 1 Introduction

Thank you very much for choosing an EIZO Monochrome Monitor.

1-1. Features

- Resolution: 5M pixels (Portrait : 2048 × 2560 dots (H × V))
- 15 Mega Sub-Pixels Display
The GS521 comes with a Sub-Pixel Drive (SPD) function capable of displaying 2048 × 7680 super high resolution by controlling each sub-pixel separately. With a 15 mega sub-pixels (MsP), sufficient information volume of the original medical image will be displayed, making the monitor suitable for rendering the outline of a detailed digital mammography image more sharply.
*15MsP viewer software needed for 15MsP display.
*10-bit simultaneous display supported when displaying 15MsP.
- Applicable to DisplayPort (applicable to 8 bit or 10 bit, not applicable to audio signals)
- 2 signal input terminals (DVI-D x 1, DisplayPort x 1)
- CAL Switch function for selecting an optimal calibration mode
“2-2. Selecting Display Mode” (page 11)
- DICOM mode (CAL Switch function) (page 29)
- The quality control software “RadiCS LE” (for Windows) used to calibrate the monitor is included
“1-3. Utility Disk” (page 7)
- The software “ScreenManager Pro for Medical” (for Windows) to adjust the screen using the mouse and keyboard is included
“1-3. Utility Disk” (page 7)
- Frame Synchronous mode supported (24.5 - 25.5Hz, 49 - 51Hz)
- HDCP (High-bandwidth Digital Content Protection)

1-2. Switches and Indicators

1. **S** switch
2. **M** switch
3. **□** switch
4. **△▽** switch
5. **○** switch
6. **⏻** switch
7. Power Indicator

Indicator status	Operation status
Green	The screen is displayed
Orange	Power saving
Off	Main Power Switch / Power Switch off

1-3. Utility Disk

An “EIZO LCD Utility Disk” (CD-ROM) is supplied with the monitor. The following table shows the disk contents and the overview of the software programs.

- **Disk contents and software overview**

The disk includes application software programs for adjustment, and User’s Manual. Refer to “Readme.txt” file on the disk for software startup procedures or file access procedures.

Item	Overview
A “Readme.txt” file	
RadiCS LE (for Windows)	RadiCS LE is quality control software used to calibrate the monitor and manage the calibration history. (A PC must be connected to the monitor with the supplied USB cable.) Refer to the description later.
RadiCS Monitor File (for Windows)	This is the monitor setting file required when using RadiCS. For information about installing this file, see the RadiCS User’s Manual.
ScreenManager Pro for Medical (for Windows)	A utility software program to control monitor adjustments from a PC using its mouse and keyboard. (A PC must be connected to the monitor with the supplied USB cable.) Refer to the description later.
User’s Manual (PDF file)	

- **To RadiCS LE or ScreenManager Pro for Medical**

Refer to the corresponding User’s Manual on the CD-ROM disk in order to install and use the software. When using this software, you will need to connect a PC to the monitor with the supplied USB cable. For more information refer to the “6-3. Making Use of USB (Universal Serial Bus)” ([page 23](#)).


1- 4. Basic Operation and Functions

Basic operation of Adjustment menu

1 Displaying Adjustment Menu


Press **O**.

The adjustment menu appears.

2 Adjusting/Setting

1. Choose a menu to adjust/set with **Δ**/**∇**, and press **O**.

2. Choose an item to adjust/set with **Δ**/**∇**, and press **O**.

3. Adjust/set the selected item with **Δ**/**∇**, and press **O**.

3 Exiting

Press **□** a few times.

The adjustment menu finishes.

Functions

The following table shows all the Adjustment menu's adjustment and setting menus.

Main Menu	Item	Adjusting/Setting
Mode 	 Brightness	"2-3. Adjusting Brightness" (page 11)
	 Mode Reset	"3-8. Restoring the Default Setting" (page 16)
PowerManager 	 Power Save	"4-1. Setting the Power Saving" (page 17)
	 Power Indicator	"4-2. Setting Power Indicator" (page 17)
Menu Settings 	 Language	"3-6. Setting Language" (page 15)
	 Menu Position	"3-7. Setting the Display Position of the Adjustment Menu" (page 15)
Tools 	 Input Selection	"6-2. Connecting Two PCs to the Monitor" (page 21)
	 Mode Preset	"3-3. Enabling/Disabling Mode Selection" (page 14)
	 Signal Info	"6-4. Displaying Monitor Information" (page 24)
	 Monitor Info	
	 All Reset	"3-8. Restoring the Default Setting" (page 16)

Chapter 2 Adjusting Screen

2-1. Setting Screen Resolution

Compatible Resolutions/Frequencies

For details on compatible resolutions, refer to “Compatible Resolutions/Frequencies” in the Setup Manual.

Setting Resolution

When you connect the monitor to the PC and find that the resolution is improper, or when you want to change the resolution, follow the procedure below.

● Windows 7

1. Right-click the mouse anywhere on the desktop except for icons.
2. From the displayed menu, click “Screen resolution”.
3. On the “Screen Resolution” dialog, select the monitor .
4. Click “Resolution” to select the desired resolution.
5. Click the [OK] button.
6. When a confirmation dialog is displayed, click [Keep changes].

● Windows Vista

1. Right-click the mouse anywhere on the desktop except for icons.
2. From the displayed menu, click “Personalize”.
3. On the “Personalization” window, click “Display Settings”.
4. On the “Display Settings” dialog, select the “Monitor” tab and select desired resolution in the “Resolution” field.
5. Click the [OK] button.
6. When a confirmation dialog is displayed, click [Yes].

● Windows XP

1. Right-click the mouse anywhere on the desktop except for icons.
2. From the displayed menu, click “Properties”.
3. When the “Display Properties” dialog is displayed, click the “Settings” tab and select desired resolution for “Screen resolution” under “Display”.
4. Click the [OK] button to close the dialog.

● Mac OS X


1. Select “System Preferences” from the Apple menu.
2. When the “System Preferences” dialog is displayed, click “Displays” for “Hardware”.
3. On the displayed dialog, select the “Display” tab and select desired resolution in the “Resolutions” field.
4. Your selection will be reflected immediately. When you are satisfied with the selected resolution, close the window.

2-2. Selecting Display Mode

CAL Switch allows you to select easily the adequate mode suitable for the monitor's application.


CAL Switch mode

Mode	Purpose
1-DICOM	Available to display in DICOM mode.
2-Native	Available for displaying images using the native characteristics of the monitor panel.
3-CAL	Displays the screen adjusted by calibration software.

Procedure

1. Press **M** .
Mode menu appears at the lower left of the screen.
2. The mode among the list is highlighted in turn each time **M** is pressed.
You can switch the current mode with Δ or ∇ while the mode menu is being displayed.

NOTE

- The Adjustment menu and the Mode menu cannot be displayed at the same time.
- You can disable a specific mode to be selected. For more information, refer to “3-3 Enabling/Disabling Mode Selection” ([page 14](#)).

2-3. Adjusting Brightness [Brightness]

The screen brightness is adjusted by changing the brightness of the backlight (Light source from the LCD back panel).

Adjustable range: 0 to 100%

Procedure

1. Choose <Mode> from the Adjustment menu, and press \circ .
2. Choose <Brightness> from the <Mode> menu, and press \circ .
3. Adjust the brightness with Δ or ∇ .
4. Press \circ to exit the adjustment.

NOTE

- You can also adjust the brightness using Δ and ∇ . (When the Adjustment menu does not appear.)

Chapter 3 Setting Monitor

3-1. Setting Orientation [Orientation]


This function changes the display orientation of the image based on the operating environment..

Attention

- Refer to the manual of the graphics board.
- The operating environments described below can be set regardless of whether or not the computer is running.

Procedure

1. Press ⏻ to turn off the monitor.
2. Press **M** holding ⏻ down for at least 2 seconds to turn on the monitor.
The Optional Settings menu appears.
3. Choose <Orientation> from the <Optional Settings> menu, and press ○ .
4. Select “Landscape”, “Portrait (SW)”, or “Portrait (HW)” with Δ or ∇ , and press ○ .

Monitor Orientation	Display Example
Landscape Select this option when using the Landscape orientation.	
Portrait (SW) Select this option when using the Portrait orientation. Graphics board utility software is used to rotate the display image 90 degrees.	
Portrait (HW) Select this option when using the Portrait orientation. The monitor function is used to rotate the display image 90 degrees.	

5. Select “Finish” with Δ or ∇ .
6. Press ○ to exit.
The Optional Settings menu is closed.

3-2. Setting the Input Signals [Signal Selection]

Signals can be selected based on your operating environment..

Attention

- A supported graphics board is required to display.
- Refer to the manual of the graphics board.

● DVI Connection (Dual Link/Single Link)

Procedure

1. Press ⏻ to turn off the monitor.
2. Press **M** holding ⏻ down for at least 2 seconds to turn on the monitor.
The Optional Settings menu appears.
3. Choose <Signal Selection> from the <Optional Settings> menu, and press ○ .
4. Select “DVI” with Δ or ∇ , and press ○ .
5. Select “Signal Format” with Δ or ∇ , and press ○ .
6. Select “Dual Link” or “Single Link” with Δ or ∇ , and press ○ .
7. Select “Sub Pixel Drive” with Δ or ∇ , and press ○ .
8. Select “On” or “Off” with Δ or ∇ , and press ○ .
9. Press \square twice, and use Δ or ∇ to select "Finish" from the <Optional Settings> menu.
10. Press ○ to exit.
The Optional Settings menu is closed.

● DVI Connection (SDG/Packed Pixel)

Procedure

1. Press ⏻ to turn off the monitor.
2. Press **M** holding ⏻ down for at least 2 seconds to turn on the monitor.
The Optional Settings menu appears.
3. Choose <Signal Selection> from the <Optional Settings> menu, and press ○ .
4. Select “DVI” with Δ or ∇ , and press ○ .
5. Select “Signal Format” with Δ or ∇ , and press ○ .
6. Select “SDG” or “Packed Pixel” with Δ or ∇ , and press ○ .
7. Select “Bit” with Δ or ∇ , and press ○ .
8. Select “8bit” or “10bit”(SDG)/“8/10bit” or “12bit”(Packed Pixel) with Δ or ∇ , and press ○ .
9. Select “Sub Pixel Drive”(SDG : 8bit) with Δ or ∇ , and press ○ .
10. Select “On” or “Off”(SDG : 8bit) with Δ or ∇ , and press ○ .
11. Press \square twice, and use Δ or ∇ to select "Finish" from the <Optional Settings> menu.
12. Press ○ to exit.
The Optional Settings menu is closed.

● DisplayPort Connection

Procedure

1. Press ⏻ to turn off the monitor.
2. Press **M** holding ⏻ down for at least 2 seconds to turn on the monitor.
The Optional Settings menu appears.
3. Choose <Signal Selection> from the <Optional Settings> menu, and press ○ .
4. Select “DisplayPort” with Δ or ∇ , and press ○ .
5. Select “Sub Pixel Drive”, and press ○ .
6. Select “On” or “Off” with Δ or ∇ , and press ○ .
7. Press \square twice, and use Δ or ∇ to select "Finish" from the <Optional Settings> menu.
8. Press ○ to exit.
The Optional Settings menu is closed.

3-3. Enabling/Disabling Mode Selection [Mode Preset]

Allows you to select the specified modes only.

Use this function when all the display modes are not available or when keeping the display mode unchanged.

Procedure

1. Choose <Tools> from the Adjustment menu, and press ○ .
2. Choose <Mode Preset> from the <Tools> menu, and press ○ .
3. Select the mode to change its settings with Δ or ∇ , and press ○ .
4. Select “On” or “Off” with Δ or ∇ .
5. Press ○ to exit the adjustment.

Attention

- You cannot set all the modes disable. Set at least one mode to “On”.
-

3-4. Locking Switches [Key Lock]

This function allows you to lock to prevent changing the adjusted/set status.

Procedure

1. Press ⏻ to turn off the monitor.
2. Press **M** holding ⏻ down for at least 2 seconds to turn on the monitor.
The Optional Settings menu appears.
3. Choose <Key Lock> from the <Optional Settings> menu, and press ○ .
4. Select “Off”, “Menu”, or “All” with Δ or ∇ , and press ○ .

Settings	Switches that can be locked
Off (Initial settings)	None (All switches are enabled)
Menu	○ switch Δ ∇ switch
All	All switches excluding ⏻

5. Select “Finish” with Δ or ∇ .
6. Press ○ to exit.
The Optional Settings menu is closed.

3-5. Setting the EIZO Logo Display [Logo]

This function allows you to display, or not, the EIZO logo.

Procedure

1. Press \odot to turn off the monitor.
2. Press **M** holding \odot down at least 2 seconds to turn on the monitor.

The Optional Settings menu appears.

3. Choose <Logo> from the <Optional Settings> menu, and press \odot .
4. Select “On” or “Off” with Δ or ∇ , and press \odot .
5. Select “Finish” with Δ or ∇ .
6. Press \odot to exit the Optional Settings.

The Optional Settings menu is closed.

3-6. Setting Language [Language]

This function allows you to select a language for the adjustment menu or displaying message.

Selectable languages

English/German/French/Spanish/Italian/Swedish/Japanese /Simplified Chinese/Traditional Chinese

Procedure

1. Choose <Menu Settings> menu from the Adjustment menu, and press \odot .
2. Choose <Language> from the <Menu Settings> menu, and press \odot .
3. Choose a language with Δ or ∇ .
4. Press \odot to exit the adjustment.

3-7. Setting the Display Position of the Adjustment Menu [Menu Position]

Adjust the menu position using the following procedure.

Procedure

1. Choose <Menu Settings> from the Adjustment menu, and press \odot .
2. Choose <Menu Position> from the <Menu Settings> menu, and press \odot .
3. Select a menu position with Δ or ∇ .
4. Press \odot to exit the adjustment.

3-8. Restoring the Default Setting

There are two types of Reset. One is used to reset the brightness adjustment only to the default settings, and the other is used to reset all the settings, except for the optional settings, to the default settings.

Attention

- After resetting, you cannot undo the operation.
-

NOTE

- For main default settings, refer to “Main default settings (factory settings)” on [page 26](#).
-

● To reset brightness adjustment values [Mode Reset]

Only the adjustment values in the current mode will revert to the default settings (factory settings).

Procedure

1. Choose <Mode> from the Adjustment menu, and press \circ .
2. Choose <Mode Reset> from the <Mode> menu and press \circ .
3. Select “Execute” with \triangle or ∇ .
4. Press \circ to exit the adjustment.

The adjustment values revert to the default settings.

● To reset all adjustments to the factory default settings [All Reset]

Reset all adjustments except for the optional settings to the factory default settings.

NOTE

- <Logo>, <Key Lock>, <Signal Selection>, and <Orientation> in the <Optional Sttings> menu are not returned to the default settings.
-

Procedure

1. Choose <Tools> from the Adjustment menu, and press \circ .
2. Choose <All Reset> from the <Tools> menu, and press \circ .
3. Select “Execute” with \triangle or ∇ .
4. Press \circ to exit the adjustment.

All settings except for the optional settings revert to the default settings.

Chapter 4 Power Saving Function

4-1. Setting the Power Saving [Power Save]

This function allows you to set the monitor into the power saving mode according to the PC status. When the monitor enters the power saving mode, no image is displayed on the screen.

Attention

- Turning off the main power switch or unplugging the power cord completely shuts off power supply to the monitor.
- Devices connected to the USB port (upstream and downstream) work when the monitor is in power saving mode or when the power switch of the monitor is Off. Therefore, power consumption of the monitor varies with connected devices even in the power saving mode.

Power Save is compliant with the following standards for the respective signal inputs.

signal inputs		Standard
Digital signal	DVI	DVI DMPM
	DisplayPort	DisplayPort Standard V1.1a

Procedure

1. Choose <PowerManager> from the Adjustment menu, and press \odot .
2. Choose <Power Save> from the <PowerManager> menu, and press \odot .
3. Select “On” or “Off” with Δ or ∇ .
4. Press \odot to exit the adjustment.

Power Saving System

The monitor enters the power saving mode in conjunction with the PC settings.

PC		Monitor Power	Indicator
Operating		Operating	Green
Power saving	STAND-BY SUSPENDED OFF	Power saving	Orange

4-2. Setting Power Indicator [Power Indicator]

The brightness of the power indicator (blue) when the screen is displayed can be adjusted (default setting is set to light up when power is turned on, and brightness is set to 4).



Procedure

1. Choose <PowerManager> from the Adjustment menu, and press \odot .
2. Choose <Power Indicator> from the <PowerManager> menu, and press \odot .
3. Select the indicator brightness “Off” or in the range from 1 to 7 with Δ or ∇ as desired.
4. Press \odot to exit the adjustment.

Chapter 5 Troubleshooting

If a problem still remains after applying the suggested remedies, contact your local dealer.

- No-picture problems → See No.1 - No.2.
- Imaging problems → See No.3 - No.7.
- Other problems → See No.8 - No.10.

Problems	Possible cause and remedy
<p>1. No picture</p> <ul style="list-style-type: none"> • Power indicator does not light. 	<ul style="list-style-type: none"> • Check whether the power cord is connected correctly. • Turn off the main power, and then turn it on again a few minutes later. • Turn the main power switch on. • Press .
<ul style="list-style-type: none"> • Power indicator lights green. • Power indicator lights orange. 	<ul style="list-style-type: none"> • Switch the input signal with S. • Operate the mouse or keyboard. • Check whether the PC is turned on.
<p>2. The message below appears.</p> <ul style="list-style-type: none"> • This message appears when no signal is input. Example: 	<p>This message appears when the signal is not input correctly even when the monitor functions properly.</p> <ul style="list-style-type: none"> • The message shown left may appear, because some PCs do not output the signal soon after power-on. • Check whether the PC is turned on. • Check whether the signal cable is connected properly. • Switch the input signal with S.
<ul style="list-style-type: none"> • The message shows that the input signal is out of the specified frequency range. Example: 	<ul style="list-style-type: none"> • Check whether the signal setting of your PC matches the resolution and the vertical frequency settings for the monitor. (refer to the Setup Manual "Setting Screen Resolution") • Reboot the PC. • Select an appropriate display mode using the graphics board's utility. Refer to the manual of the graphics board for details. fD : Dot Clock fH : Horizontal Frequency fV : Vertical Frequency
<p>3. The screen is too bright or too dark.</p>	<ul style="list-style-type: none"> • Adjust using <Brightness>. (The LCD monitor backlight has a fixed life span. When the screen becomes dark or begins to flicker, contact your local dealer.)
<p>4. Characters are blurred</p>	<ul style="list-style-type: none"> • Check whether the signal setting of your PC matches the resolution and the vertical frequency settings for the monitor. (refer to the Setup Manual "Setting Screen Resolution")
<p>5. Afterimages appear</p>	<ul style="list-style-type: none"> • Afterimages are particular to LCD monitors. Avoid displaying the same image for a long time. • Use the screen saver or power save function to avoid displaying the same image for extended periods of time.

Problems	Possible cause and remedy
6. Interference patterns or pressure marks remain on the screen.	<ul style="list-style-type: none"> • Leave the monitor with a white or black screen. The symptom may disappear.
7. Noise appears on the screen.	<ul style="list-style-type: none"> • When entering the signals of HDCP system, the normal images may not be displayed immediately.
8. The adjustment menu does not appear.	<ul style="list-style-type: none"> • Check whether the operation lock function works. (page 14)
9. Mode menu does not appear.	<ul style="list-style-type: none"> • Check whether the operation lock function works. (page 14)
10. The monitor connected with the USB cable is not detected. / USB devices connected to the monitor does not work.	<ul style="list-style-type: none"> • Check whether the USB cable is connected correctly. (page 23) • Change the USB port to another one. If the PC or peripheral devices works correctly by changing the USB port, contact your local dealer. (Refer to the manual of the PC for details.) • Reboot the PC. • If the peripheral devices work correctly when the PC and peripheral devices are connected directly, please contact your local dealer. • Check whether the PC and OS are USB compliant. (For USB compliance of the respective devices, consult their manufacturers.) • Check the PC's BIOS setting for USB when using Windows. (Refer to the manual of the PC for details.)

Chapter 6 Reference

6-1. Attaching an Arm

The stand can be removed and replaced with an arm (or another stand) to be attached to the monitor. Use an arm or stand of EIZO option.

Attention

- When attaching an arm or stand, follow the instructions of their user's manual.
- When using another manufacturer's arm or stand, confirm the following in advance and select one conforming to the VESA standard.
 - Clearance between the screw holes: 100 mm × 100 mm
 - Thickness of plate: 2.6 mm
 - Strong enough to support weight of the monitor unit (except the stand) and attachments such as cables.
- When using an arm or stand, attach it to meet the following tilt angles of the monitor.
 - Up 45 degrees, down 45 degrees (horizontal display, and vertical display rotated 90 degrees counter-clockwise)
- Connect the cables after attaching an arm or a stand
- Do not adjust the height of the stand removed from the monitor. When adjusted without attached monitor, it causes the personnel injury or the damaged to the stand.
- Since the monitor and arm or stand are so heavy, dropping them may result in injury or equipment damage.

[Attachment Procedure]


1 Lay the LCD monitor on a soft cloth spread over on a stable surface with the panel surface facing down.

2 Remove the stand.

Prepare a screwdriver. Unscrew the four screws securing the unit and the stand with the screwdriver.

3 Attach the monitor to the arm or stand.


Secure the monitor to the arm or stand using the screws specified in the user's manual of the arm or stand.

6-2. Connecting Two PCs to the Monitor

Two PCs can be connected to the monitor through the DVI-D and the DisplayPort connector on the back of the monitor.


Connection examples

● To switch the input signal

The input signal switches each time **S** is pressed.

When the signal is switched, the active input port name appears at the top right corner of the screen.

● To set input signal selection [Input Selection]

Setting	Function
Auto	The monitor recognizes the connector through which PC signals are input. When a PC is turned off or enters the powersaving mode, the monitor automatically displays another signal.
Manual	The monitor detects only the PC's signals currently displaying automatically. Select an active input signal with S .

Procedure

1. Choose <Tools> from the Adjustment menu, and press **O** .
2. Choose <Input Selection> from the <Tools> menu, and press **O** .
The <Input Selection> menu appears.
3. Select “Auto” or “Manual” with **Δ** or **∇** .
4. Press **O** to exit the adjustment.

Attention

- When “Auto” is selected for <Input Signal>, the monitor’s power-saving function works only when all PCs are in the power-saving mode.
-

6-3. Making Use of USB (Universal Serial Bus)

This monitor has a hub compatible with USB. Connected to a PC compatible with USB or another USB hub, this monitor functions as a USB hub allowing connection to peripheral USB devices.

● Required System Environment


1. A PC equipped with a USB port or another USB hub connected to a USB compatible PC
2. Windows 2000/XP/Vista/7 or Mac OS 9.2.2 and Mac OS X 10.2 or later
3. EIZO USB cable (MD-C93)

Attention

- This monitor may not work depending on PC, OS or peripheral devices to be used. For USB compatibility of peripheral devices, contact their manufactures.
- Devices connected to the USB port (upstream and downstream) work when the monitor is in power saving mode or when the power switch of the monitor is Off. Therefore, power consumption of the monitor varies with connected devices even in the power saving mode.
- When the main power switch is Off, device connected to the USB port will not operate.

● Connection Procedure (Setup of USB Function)

1. Connect the monitor first to a PC using the signal cable, and run the PC.
2. Connect the supplied USB cable between the downstream USB port of a USB compatible PC (or another USB hub) and the monitor's upstream USB port.
The USB hub function is set up automatically upon connection of the USB cable.
3. Connect the peripheral USB device to the USB port (downstream) of the monitor.


6-4. Displaying Monitor Information

● Displaying signal information [Signal Info]

This function displays the information about the current input signals displayed.

1. Choose <Tools> from the Adjustment menu, and press \odot .
 2. Choose <Signal Info> from the <Tools> menu, and press \odot .
- The Signal Info screen appears.

(Example)


● Displaying monitor information [Monitor Info]

This function displays the information about the monitor.

1. Choose <Tools> from the Adjustment menu, and press \odot .
 2. Choose <Monitor Info> from the <Tools> menu, and press \odot .
- The Monitor Info screen appears.

(Example)

Attention

- The usage time is not always "0" due to factory inspection when you purchase the monitor.
-

6-5. Specifications

LCD Panel	Size	21.3-inch (540 mm) TFT Monochrome LDC panel
	Surface treatment (GS521-CL/BL)	Anti-Glare
	Surface treatment (GS521-CLG)	Glare
	Surface hardness	3H
	Viewing angle	Viewing angle : Horizontal 170°, Vertical 170° (CR:10 or more)
	Dot Pitch	0.165 mm
	Response Time	Approx. 50 ms
Horizontal Scan Frequency		31-135 kHz
Vertical Scan Frequency		19 ~ 61Hz (VGA TEXT : 69 ~ 71Hz / QSXGA : 19 ~ 51Hz)
Resolution		5M pixels (Portrait: 2048 × 2560 dots (H × V)) 15M Sub-Pixels (Portrait: 2048 × 7680 dots (H × V) ([Sub Pixel Drive] : "On"))
Max. Dot Clock		290 MHz
Gray Scale		4096 steps of 13771 steps (12 bit signal input, [Sub Pixel Drive] : "Off") 1024 steps of 3571 steps (10 bit signal input, [Sub Pixel Drive] : "On")
Recommended Brightness		500 cd/m ² (approx.65%)
Display Area (H × V)		337.9 mm × 422.4 mm
Power Supply		100-120 VAC ±10%, 50/60Hz 1.1 – 0.9 A 200-240 VAC ±10%, 50/60Hz 0.6 – 0.5 A
Power Consumption	Screen Display On	100W or less (with USB load) 90W or less (without USB load)
	Power saving mode	1.0W or less (for DVI-D single signal input, without USB load, [Input Signal] : "Manual")
	Power switch Off	0.8W or less (without USB load)
	Main Power switch Off	0W
Input Signal Connector		DVI-D connector (Applicable to HDCP) × 1
		DisplayPort (Standard V1.1a, applicable to HDCP)
Digital Signal (DVI) Transmission System		TMDS (Single Link/Dual Link)
Plug & Play		Digital (DVI-D) : VESA DDC 2B / EDID structure 1.3 Digital (DisplayPort) : VESA DisplayPort / EDID structure 1.4
Dimensions	Main unit	388.0 mm (15.3 inch) × 512.5-594.5 mm (20.2 – 23.4 inch) × 208.5 mm (8.2 inch)
	Main unit (without stand)	388.0 mm (15.3 inch) × 480 mm (18.9 inch) × 88 mm (3.46 inch)
Mass	Main unit	Approx. 10.1 kg (22.3 lbs.)
	Main unit (without stand)	Approx.7.1 kg (15.7 lbs.)
Movable range	Height adjustable stand	Tilt: 40° Up, 0° Down Swivel: 35° Right, 35° Left Adjustable height: 82 mm (3.23 inch) Rotation: 90° (counter-clockwise)

Environmental Conditions	Temperature	Operating: 0 °C - 35 °C (32 °F - 95 °F) Transportation/Storage: -20 °C - 60 °C (-4 °F - 140 °F)
	Humidity	Operating: 20% - 80% R.H. (no condensation) Transportation/Storage: 10% - 80% R.H. (no condensation)
	Pressure	Operating: 700 to 1,060 hPa Transportation/Storage: 200 to 1,060 hPa
USB	Standard	USB Specification Revision 2.0
	Port	Upstream port × 1, Downstream port × 2
	Communication Speed	480 Mbps (high), 12 Mbps (full), 1.5 Mbps (low)
	Supply current	Downstream: Max. 500mA/1 port
Classification of Equipment		Type of protection against electric shock : Class I EMC class: EN60601-1-2:2007 Group 1 Class B Classification of medical device (MDD 93/42/EEC) : Class I Mode of operation: Continuous IP Class: IPX0


Main default settings (factory settings)

CAL Switch	1-DICOM
Power Save	On
Input Selection	Auto
Logo	Off
Menu Position	Center
Language	English
Signal Selection*	Single Link
Orientation*	Portrait (HW)

* These functions cannot be initialized with the All Reset function ([page 16](#)). For information about setting these functions, refer to "3-2. Setting the Input Signals" ([page 13](#)) and "3-1. Setting Orientation" ([page 12](#)).


Outside Dimensions

Unit: mm (inch)

Connector Pin Assignment


- DVI-D connector

Pin No.	Signal	Pin No.	Signal	Pin No.	Signal
1	T.M.D.S. Data 2-	9	T.M.D.S. Data1-	17	T.M.D.S. Data0-
2	T.M.D.S. Data 2+	10	T.M.D.S. Data1+	18	T.M.D.S. Data0+
3	T.M.D.S. Data2/4 Shield	11	T.M.D.S. Data1/3 Shield	19	T.M.D.S. Data0/5 Shield
4	T.M.D.S. Data 4-	12	T.M.D.S. Data 3-	20	T.M.D.S. Data 5-
5	T.M.D.S. Data 4+	13	T.M.D.S. Data 3+	21	T.M.D.S. Data 5+
6	DDC Clock (SCL)	14	+5V Power	22	T.M.D.S. Clock shield
7	DDC Data (SDA)	15	Ground (return for +5V, Hsync and Vsync)	23	T.M.D.S. Clock+
8	NC*	16	Hot Plug Detect	24	T.M.D.S. Clock-


(NC*: No Connection)

- DisplayPort connector

Pin No.	Signal	Pin No.	Signal	Pin No.	Signal
1	ML Lane3-	8	Ground	15	AUX CH+
2	Ground	9	ML Lane1+	16	Ground
3	ML Lane3+	10	ML Lane0-	17	AUX CH-
4	ML Lane2-	11	Ground	18	Hot Plug Detect
5	Ground	12	ML Lane0+	19	Return
6	ML Lane2+	13	CONFIG1	20	DP PWR
7	ML Lane1-	14	CONFIG2		

- USB port

Contact No.	Signal	Remarks
1	VCC	Cable power
2	- Data	Serial data
3	+ Data	Serial data
4	Ground	Cable ground

Accessories List

Calibration Kit	EIZO "RadiCS UX1" Ver.3.4.0 or later
Network QC Management Software	EIZO "RadiNET Pro" Ver.3.4.0 or later
Cleaning Kit	EIZO "ScreenCleaner"

For the latest information about the accessories and information about the latest compatible graphics board, refer to our web site.

<http://www.eizo.com>

6-6. Glossary

DICOM (Digital Imaging and Communication in Medicine)

DICOM is a standard for digital imaging and communication for medical use developed by American College of Radiology and National Electric Manufacturers Association.

DisplayPort

DisplayPort is a next-generation digital AV interface that allows connection of the PC, audio, imaging devices, etc. to the monitor. One cable can transfer sound with images.

DVI (Digital Visual Interface)

DVI is a digital interface standard. DVI allows direct transmission of the PC's digital data without loss. This adopts the TMDS transmission system and DVI connectors. There are two types of DVI connectors. One is a DVI-D connector for digital signal input only. The other is a DVI-I connector for both digital and analog signal inputs.

DVI DMPM (DVI Digital Monitor Power Management)

DVI DMPM is a digital interface power-saving function. The "Monitor ON (operating mode)" and "Active Off (power-saving mode)" are indispensable for DVI DMPM as the monitor's power mode.

HDCP (High-bandwidth Digital Content Protection)

Digital signal coding system developed to copy-protect the digital contents, such as video, music, etc. This helps to transmit the digital contents safely by coding the digital contents sent via DVI terminal on the output side and decoding them on the input side. Any digital contents cannot be reproduced if both of the equipments on the output and input sides are not applicable to HDCP system.

Resolution

The LCD panel consists of numerous pixels of specified size, which are illuminated to form images. This monitor consists of 2048 horizontal pixels and 2560 vertical pixels. At a resolution of 2048 x 2560, all pixels are illuminated as a full screen (1:1).

TMDS (Transition Minimized Differential Signaling)

A signal transmission system for digital interface.

For U.S.A. , Canada, etc. (rated 100-120 Vac) Only

FCC Declaration of Conformity

We, the Responsible Party

EIZO NANAOTECHNOLOGIES INC.

5710 Warland Drive, Cypress, CA 90630

Phone: (562) 431-5011

declare that the product

Trade name: EIZO

Model: RadiForce GS521

is in conformity with Part 15 of the FCC Rules. Operation of this product is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures.

- * Reorient or relocate the receiving antenna.
- * Increase the separation between the equipment and receiver.
- * Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- * Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note

Use the attached specified cable below or EIZO signal cable with this monitor so as to keep interference within the limits of a Class B digital device.

- AC Cord
- Shielded Signal Cable (Enclosed)

Canadian Notice

This Class B digital apparatus complies with Canadian ICES-003.

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EMC Information

Essential performance of RadiForce series is to display images and operate functions normally.

CAUTION

The RadiForce series requires special precautions regarding EMC and need to be installed, put into service and used according to the following information.

Do not use any cables other than the cables that provided or specified by us.

Using other cables may cause the increase of emission or decrease of immunity.

Do not put any portable and mobile RF communications equipment close to the RadiForce series. Doing so may affect the RadiForce series.

The RadiForce series should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

Anyone who connects additional equipment to the signal input part or signal output parts, configuring a medical system, is responsible that the system complies with the requirements of IEC/EN60601-1-2.

Guidance and manufacturer's declaration - electromagnetic emissions

The RadiForce series is intended for use in the electromagnetic environment specified below.

The customer or the user of the RadiForce series should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11/EN55011	Group 1	The RadiForce series uses RF energy only for its internal function. Therefore, its RF emission are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11/EN55011	Class B	The RadiForce series is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC/EN61000-3-2	Class D	
Voltage fluctuations / flicker emissions IEC/EN61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity

The RadiForce series is intended for use in the electromagnetic environment specified below.

The customer or the user of the RadiForce series should assure that it is used in such an environment.

Immunity test	IEC/EN60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC/EN61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC/EN61000-4-5	±1kV line(s) to line(s) ±2kV line(s) to earth	±1kV line(s) to line(s) ±2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC/EN61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the RadiForce series requires continued operation during power mains interruptions, it is recommended that the RadiForce series be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC/EN61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity			
The RadiForce series is intended for use in the electromagnetic environment specified below. The customer or the user of the RadiForce series should assure that it is used in such an environment.			
Immunity test	IEC/EN60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC/EN61000-4-6	3Vrms 150kHz to 80MHz	3Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the RadiForce series, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended Separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$, 80MHz to 800MHz $d = 2.3 \sqrt{P}$, 800MHz to 2.5GHz Where "P" is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and "d" is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol. 
Radiated RF IEC/EN61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the RadiForce series is used exceeds the applicable RF compliance level above, the RadiForce series should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the RadiForce series.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the RadiForce Series

The RadiForce series is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the RadiForce series can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the RadiForce series as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = 1.2 \sqrt{P}$	80MHz to 800MHz $d = 1.2 \sqrt{P}$	800MHz to 2.5GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance “d” in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where “P” is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Cable length		
Power Cord :	Accessory	2.0m
Signal Cable (DD200-DL) :	Accessory	2.0m
Signal Cable (PP200) :	Accessory	2.0m
USB Cable (MD-C93) :	Accessory	1.8m

LIMITED WARRANTY

EIZO NANO CORPORATION (hereinafter referred to as “EIZO”) and distributors authorized by EIZO (hereinafter referred to as the “Distributors”) warrant, subject to and in accordance with the terms of this limited warranty (hereinafter referred to as the “Warranty”), to the original purchaser (hereinafter referred to as the “Original Purchaser”) who purchased the product specified in this document (hereinafter referred to as the “Product”) from EIZO or Distributors, that EIZO and Distributors shall, at their sole discretion, either repair or replace the Product at no charge if the Original Purchaser becomes aware within the Warranty Period (defined below) that (i) the Product malfunctions or is damaged in the course of normal use of the Product in accordance with the description in the instruction manual attached to the Product (hereinafter referred to as the “User’s Manual”), or (ii) the LCD panel and brightness of the Product cannot maintain the recommended brightness specified in the User’s Manual in the course of normal use of the Product in accordance with the description of the User’s Manual.

The period of this Warranty is five (5) years from the date of purchase of the Product (hereinafter referred to as the “Warranty Period”), provided that the Warranty Period for the LCD panel is limited to five (5) years from the date of purchase of the Product subject to the usage time being less than or equal to 30,000 hours. The brightness of the Product shall be warranted only if the Product has been used within the recommended brightness of described in the User’s Manual and the Warranty Period for the brightness is limited to five (5) years from the date of purchase of the Product subject to the usage time being less than or equal to 30,000 hours. EIZO and Distributors shall bear no liability or obligation with regard to the Product in relation to the Original Purchaser or any third parties other than as provided under this Warranty.

EIZO and Distributors will cease to hold or store any parts of the Product upon expiration of seven (7) years after the production of such parts is discontinued. In repairing the monitor, EIZO and Distributors will use renewal parts which comply with our QC standards.

The Warranty is valid only in the countries or territories where the Distributors are located. The Warranty does not restrict any legal rights of the Original Purchaser.

Notwithstanding any other provision of this Warranty, EIZO and Distributors shall have no obligation under this Warranty whatsoever in any of the cases as set forth below:

- (a) Any defect of the Product caused by freight damage, modification, alteration, abuse, misuse, accident, incorrect installation, disaster, faulty maintenance and/or improper repair by third party other than EIZO and Distributors;
- (b) Any incompatibility of the Product due to possible technical innovations and/or regulations;
- (c) Any deterioration of the sensor;
- (d) Any deterioration of display performance caused by the deterioration of expendable parts such as the LCD panel and/or backlight, etc. (e.g. changes in color, changes in color uniformity, defects in pixels including burnt pixels, etc.);
- (e) Any deterioration of the Product caused by the use at higher brightness than the recommended brightness described in the User’s Manual;
- (f) Any defect of the Product caused by external equipment;
- (g) Any defect of the Product on which the original serial number has been altered or removed;
- (h) Any normal deterioration of the product, particularly that of consumables, accessories, and/or attachments (e.g. buttons, rotating parts, cables, User’s Manual, etc.); and
- (i) Any deformation, discoloration, and/or warp of the exterior of the product including that of the surface of the LCD panel.

To obtain service under the Warranty, the Original Purchaser must deliver the Product, freight prepaid, in its original package or other adequate package affording an equal degree of protection, assuming the risk of damage and/or loss in transit, to the local Distributor. The Original Purchaser must present proof of purchase of the Product and the date of such purchase when requesting services under the Warranty.

The Warranty Period for any replaced and/or repaired product under this Warranty shall expire at the end of the original Warranty Period.

EIZO OR DISTRIBUTORS ARE NOT RESPONSIBLE FOR ANY DAMAGE TO, OR LOSS OF, DATA OR OTHER INFORMATION STORED IN ANY MEDIA OR ANY PART OF ANY PRODUCT RETURNED TO EIZO OR DISTRIBUTORS FOR REPAIR.

EIZO AND DISTRIBUTORS MAKE NO FURTHER WARRANTIES, EXPRESSED OR IMPLIED, WITH RESPECT TO THE PRODUCT AND ITS QUALITY, PERFORMANCE, MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR USE. IN NO EVENT SHALL EIZO OR DISTRIBUTORS BE LIABLE FOR ANY INCIDENTAL, INDIRECT, SPECIAL, CONSEQUENTIAL OR OTHER DAMAGE WHATSOEVER (INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF PROFIT, BUSINESS INTERRUPTION, LOSS OF BUSINESS INFORMATION, OR ANY OTHER PECUNIARY LOSS) ARISING OUT OF THE USE OR INABILITY TO USE THE PRODUCT OR IN ANY CONNECTION WITH THE PRODUCT, WHETHER BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE, EVEN IF EIZO OR DISTRIBUTORS HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. THIS EXCLUSION ALSO INCLUDES ANY LIABILITY WHICH MAY ARISE OUT OF THIRD PARTY CLAIMS AGAINST THE ORIGINAL PURCHASER. THE ESSENCE OF THIS PROVISION IS TO LIMIT THE POTENTIAL LIABILITY OF EIZO AND DISTRIBUTORS ARISING OUT OF THIS LIMITED WARRANTY AND/OR SALES.

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