



## EIZO Upgrades RadiCS Quality Control Software to Comply with ACR Practice Guideline

Hakusan, Japan, November, 1, 2007 – Eizo Nanao Corporation (TSE: 6737) today announced it is upgrading its RadiCS quality control software to version 3.1.0 to comply with the American College of Radiology (ACR) “Practice Guideline for Determinants of Image Quality in Digital Mammography released in October 2007.”



This guideline was developed by the ACR in cooperation with the American Association of Physicists in Medicine (AAPM), Radiological Society of North America (RSNA), and Society for Imaging Informatics in Medicine (SIIM). The guideline is an educational tool to provide radiologists, medical physicists and other support personnel with in-depth knowledge of the issues pertinent to digital mammography image quality. In order to explain monitor quality control (QC), the guideline frequently refers to AAPM Task Group 18 documentation. Utilizing its expertise in monitor QC, EIZO is among the first to implement this guideline into its quality control software.

As with previous versions, RadiCS version 3.1.0 enables you to perform pattern, luminance, grayscale and uniformity checks complying with AAPM TG18, DIN V 6868-57, and other QC standards. In addition, it now complies with QC standards for digital mammography monitors such as PAS 1054, EUREF European Guidelines, and the ACR Practice Guideline presenting healthcare and IT professionals in the mammography field with a cost-effective and time-efficient way to ensure stable, consistent image quality from all their EIZO monitors, which ultimately leads to the improvement in the quality of medical services.

### Availability

Availability varies by country so please contact the EIZO subsidiary or distributor in your country for details. Visitors to the EIZO booth at RSNA (Radiological Society of North America) 2007 in Chicago, Illinois, USA can have a first look at the RadiCS Version 3.1.0 in the EIZO booth located in the North building, Hall B, Booth 8534. RSNA runs from November 25-29.

### About RadiCS Quality Control Software

RadiCS provides total support for the quality maintenance and control of client monitors, covering everything from calibration to acceptance and constancy tests, calibration history and asset management. Complying with AAPM, DIN, IEC and other standards, RadiCS enables precise QC with easy-to-follow procedures.

### About EIZO

Eizo Nanao Corporation is a leading global manufacturer of high-end visual display products with a wide range of LCD monitors. The image quality, long-term reliability, and innovative features of EIZO monitors make them the products of choice in many financial trading rooms, hospitals, back offices, and design studios throughout the world. EIZO is based in Japan and represented in over fifty countries by a network of exclusive distributors.



For more information, please contact:

Eizo Nanao Corporation  
153 Shimokashiwano  
Hakusan, Ishikawa 924-8566  
Japan  
Phone: +81 76 277-6792  
Fax: +81 76 277-6793  
<http://www.radiforce.com>

###

*RadiCS is a trademark, and RadiForce and EIZO are registered trademarks of Eizo Nanao Corporation.*